

---

## 1 API Overview

### Base URL

https://hospital-domain/api/v1

### Auth

- **JWT** for staff users
- **mTLS** or **device tokens** for hardware devices
- **Role-based access control (RBAC)**

### Response format

```
{  
  "success": true,  
  "data": { ... },  
  "error": null  
}
```

---

## 2 Authentication APIs

### POST /auth/login

Authenticate staff user.

### Body

```
{  
  "username": "nurse1",  
  "password": "*****"  
}
```

### Response

```
{  
  "token": "jwt-token",  
  "expires_in": 3600  
}
```

---

### Mother & Infant Management

#### POST /mothers

Register a mother.

##### Body

```
{  
  "hospital_mrn": "MRN12345",  
  "first_name": "Sara",  
  "last_name": "Ali",  
  "national_id": "1234567890",  
  "phone": "0500000000"  
}
```

---

#### POST /infants

Register an infant.

##### Body

```
{  
  "hospital_infant_no": "INF-2024-001",  
  "gender": "FEMALE",  
  "date_of_birth": "2024-01-01T10:30:00Z",  
  "mother_id": "uuid"  
}
```

---

#### GET /infants/{infantId}

Retrieve infant details.

---

### Tag Assignment & Pairing

#### POST /tags/infant/assign

Assign infant tag.

**Body**

```
{  
  "infant_id": "uuid",  
  "tag_uid": "TAG123ABC"  
}
```

---

**POST /tags/mother/assign**

Assign mother tag.

---

**POST /pairings**

Pair infant + mother + tags.

**Body**

```
{  
  "infant_id": "uuid",  
  "mother_id": "uuid",  
  "infant_tag_uid": "TAG-INF-001",  
  "mother_tag_uid": "TAG-MOM-001"  
}
```

---

**GET /pairings/{infantId}**

Retrieve active pairing.

---

**5 RTLS & Location APIs**

**GET /location/tag/{tagUid}**

Returns current zone + last seen time.

**Response**

```
{
```

```
"tag_uid": "TAG-INF-001",  
"zone_id": "uuid",  
"zone_name": "Nursery Room",  
"last_seen_at": "2024-01-01T12:00:00Z"  
}
```

---

### **GET /location/infant/{infantId}**

Returns infant's current location.

---

### **POST /rtls/readerEvent (*device* → *server*)**

Reader sends tag sighting.

#### **Body**

```
{  
  "reader_id": "uuid",  
  "tag_uid": "TAG-INF-001",  
  "rssi": -55,  
  "timestamp": "2024-01-01T12:00:00Z"  
}
```

---

## **6 Gate Movement Authorization**

### **POST /gate/authorizeMovement**

Gate terminal requests permission.

#### **Body**

```
{  
  "infant_tag_uid": "TAG-INF-001",  
  "mother_tag_uid": "TAG-MOM-001",  
  "staff_id": "uuid",  
  "gate_id": "uuid",  
}
```

```
"reason": "ROUTINE_TRANSFER"
}
```

## Response

```
{
  "authorized": true,
  "from_zone": "Nursery",
  "to_zone": "NICU",
  "movement_id": "uuid"
}
```

---

## GET /gate/movements

Filter movement logs.

Query params:

?infant\_id=&mother\_id=&staff\_id=&from\_date=&to\_date=

---

## 7 Tamper & Security Events

### POST /events/tamper (*device* → *server*)

#### Body

```
{
  "tag_uid": "TAG-INF-001",
  "device_id": "uuid",
  "timestamp": "2024-01-01T12:00:00Z"
}
```

---

## GET /events

List events.

Query params:

?type=TAMPER&severity=CRITICAL&from\_date=&to\_date=

---

## **8 Alarm Control APIs**

### **POST /alarms/raise**

Raise alarm manually or from event.

#### **Body**

```
{  
  "zone_id": "uuid",  
  "event_id": "uuid"  
}
```

---

### **POST /alarms/silence**

Silence active alarm.

---

### **POST /alarms/reset**

Reset alarm state.

---

### **GET /alarms/status**

Returns all active alarms.

---

## **9 Biometric APIs**

### **POST /biometric/enrollInfant**

Enroll footprint template.

#### **Body**

```
{  
  "infant_id": "uuid",  
  "template_format": "ISO-FMR",  
  "template_data": "base64-encoded",  
  "quality_score": 92  
}
```

```
}
```

---

### **POST /biometric/verifyInfant**

Verify infant footprint.

#### **Body**

```
{  
  "infant_id": "uuid",  
  "template_data": "base64-encoded"  
}
```

---

### **POST /biometric/checkDuplicate**

Check if footprint already exists.

---

## **10 Device Management APIs**

### **POST /devices/register**

Register hardware device.

#### **Body**

```
{  
  "device_type": "READER",  
  "serial_number": "SN12345",  
  "firmware_version": "1.0.0"  
}
```

---

### **POST /devices/heartbeat**

Device health ping.

#### **Body**

```
{  
  "device_id": "uuid",
```

```
"status": "OK",  
"battery": 85,  
"temperature": 40  
}
```

---

### **GET /devices**

List all devices.

---

## **1 1 Configuration APIs**

### **GET /config/zones**

Returns hospital topology.

---

### **POST /config/zones**

Create new zone.

---

### **POST /config/gates**

Create or update gate.

---

## **1 2 Audit APIs**

### **GET /audit**

Query audit logs.

Query params:

?user\_id=&action\_type=&entity\_type=&from\_date=&to\_date=

---

## **1 3 WebSocket Events**

### **/ws/alerts**

Push notifications for:

- Tamper events



- Unauthorized movement
- Alarms
- Infant location changes

### **/ws/dashboard**

Real-time dashboard updates.

---

#### **Error Response Format**

```
{  
  "success": false,  
  "error": {  
    "code": "INVALID_PAIRING",  
    "message": "Mother tag does not match infant pairing."  
  }  
}
```

---